PXDENTAL SA

BLANC	ALLIAGE DENTAIRE CÉRAMO MÉTALLIQUE
WEISS	DENTAL AUFBRENNLEGIERUNGEN
BIANCO	LEGA DENTALE PER CERAMICA
WHITE	DENTAL BONDING ALLOY
BLANCO	ALECION DENTAL CERAMO-METALLICA

% = Au: 54.2, Pd: 31.0, Ag: 4.8, In: 9.0, Ru, Ga: < 1%

14.6 g/cm³
1160-1260° C
14.3/14.5 µm/m² C

ISO 9693

Fabriqué par / Hergestellt von / Fabricato da /
Made by / Fabricado por

PX DENTAL SA

Champs-Montants 16a, 2074 Marin
T. +41 (0)32 924 21 20 F +41 (0)32 924 21 29
www.pxdental.com

PXDENTAL SA
CH - 2074 MARIN

Mode d'emploi

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1. Waxing up, spruing, investing

- Follow the general prescriptions for casting alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- All regular or phosphate-bonded investments for precious metal alloys may be used. For the alloy PX BONDALLOY we recommend to use the PX Extravest investment.

2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 850° C for a soaking period of 30-45 min.

3. Crucibles and casting buttons

- PX BONDALLOY may be cast in a graphite, quartz or ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

4. Melting

- For PX BONDALLOY casting temperature is 1440° C with a ceramic crucible.
- The use of a pinch of boric acid is recommended.
- Bench cool the casting ring.

5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Finishing: use abrasive grinding points of your choice. Clean the frameworks in an ultrasonic unit or with a steam cleaner.

6. Oxidation

- 10 min without vacuum at 950°C.

7. Porcelain application

- Porcelain applications are subject to the manufacturer's instructions.
- PX BONDALLOY is compatible with most known porcelain brands. (high fusing).

8. Soldering

- Pre-bond : use PX PRESOLD 1110 at 1110° C with flux.
- Post-bond : use PX POSTSOLD 800 at 850° C with flux.

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