N° 7151

PX EXTRAPALL 2

**PXDENTAL SA**

|        |                                    |
|--------|------------------------------------|
| BLANC  | ALLIAGE DENTAIRE CÉRAMO MÉTALLIQUE |
| WEISS  | DENTAL AUFBRENNLEGIERUNGEN         |
| BIANCO | LEGA DENTALE PER CERAMICA          |
| WHITE  | DENTAL BONDING ALLOY               |
| BLANCO | ALECION DENTAL CERAMO-METALLICA    |

% = Au: 2.0, Pd: 78.5, Cu: 6.9, In: 4.5, Sn: 2.0,  
Ga: 5.5, Ru, Zn: < 1%

11.4 g/cm<sup>3</sup>  
1150-1270° C  
13.9/14.2 µm/m° C

ISO 9693

Fabriqué par / Hergestellt von / Fabricato da /  
Made by / Fabricado por

**PX DENTAL SA**

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**1. Waxing up, spruing, investing**

- Follow the general prescriptions for bonding alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- Use only a high-heat phosphate-bonded investment material.

**2. Preheating**

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 850° C for a soaking period of 20-45 min.

**3. Crucibles and casting buttons**

- Use a ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

**4. Melting**

- For PX EXTRAPALL 2 casting temperature is 1450° C.
- The use of crystallized boric acid is recommended.
- Bench cool the casting ring.

**5. Treatment of the casting**

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and preceding oxidation clean the frameworks in an ultrasonic unit or with a steam cleaner.

**6. Oxidation**

- 10 min without vacuum at 950° C.

**7. Porcelain application**

- Porcelain applications are subject to the manufacturers instructions.
- PX EXTRAPALL 2 is compatible with most known porcelain brands (high fusing).

**8. Soldering**

- Pre-bond: use PX PRESOLD 1110 solder at 1110° C with a flux.
- Post-bond: use PX POSTSOLD 800 solder at 850° C with a flux.

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