N° 7301

PX BIO GOLD LF

PX DENTAL SA

JAUNE	ALLIAGE DENTAIRE UNIVERSEL
GELB	DENTAL UNIVERSALLEGIERUNGEN
GIALLO	LEGA DENTALE UNIVERSAL
YELLOW	DENTAL UNIVERSAL ALLOY
AMARILLO	ALEACIÓN DENTAL UNIVERSAL

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16.7 g/cm<sup>3</sup>

970-1060° C

15.6/15.9 µm/m° C

ISO 9693 /22674

Fabriqu  par / Hergestellt von / Fabricato da /  
Made by / Fabricado por

PX DENTAL SA

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## 1. Waxing up, spruing, investing

- Follow the general prescriptions for bonding alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- Use only a high-heat phosphate-bonded investment material. For the alloy PX BIO GOLD LF we recommend to use the PX Extravest investment.

## 2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 700° C for a soaking period of 20-45 min.

## 3. Crucibles and casting buttons

- May be cast in a graphite or ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

## 4. Melting

- For PX BIO GOLD LF casting temperature is 1210° C with a graphite crucible, 1260° C with a ceramic crucible.
- The use of crystallized boric acid is recommended.
- Bench cool the casting ring.

## 5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and preceding oxidation clean the frameworks in an ultrasonic unit or with a steam cleaner.

## 6. Oxidation

- 10 min without vacuum at 800° C.

## 7. Porcelain application

- Porcelain applications are subject to the manufacturers instructions.
- PX BIO GOLD LF is compatible with low fusing ceramics (like Duceragold). Please observe the required coefficient of expansion.

## 8. Soldering

- Pre-bond: use PX Sold LF 950 solder at 950° C with a flux.
- Post-bond: use PX Sold 700 solder at 750° C with a flux.

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