N° 7302

PX GOLD SUPRA

PX DENTAL SA

JAUNE	ALLIAGE DENTAIRE UNIVERSEL
GELB	DENTAL UNIVERSALLEGIERUNGEN
GIALLO	LEGA DENTALE UNIVERSAL
YELLOW	DENTAL UNIVERSAL ALLOY
AMARILLO	ALEACIÓN DENTAL UNIVERSAL

% = Au: 73.0, Pt: 3.0, Pd: 6.5, Ag: 14.0, Zn: 2.8, Ir, Sn: < 1%

16.1 g/cm³

1030-1080° C

16.3/16.6 µm/m° C

ISO 9693 /22674

Fabriqu  par / Hergestellt von / Fabricato da /
Made by / Fabricado por

PX DENTAL SA

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1. Waxing up, spruing, investing

- Follow the general prescriptions for bonding alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- Use only a high-heat phosphate-bonded investment material. For the alloy PX GOLD SUPRA we recommend to use the PX Extravest investment.

2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 750° C for a soaking period of 20-45 min.

3. Crucibles and casting buttons

- May be cast in a graphite or ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

4. Melting

- For PX GOLD SUPRA casting temperature is 1230° C with a graphite crucible, 1280° C with a ceramic crucible.
- The use of crystallized boric acid is recommended.
- Bench cool the casting ring.

5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and preceding oxidation clean the frameworks in an ultrasonic unit or with a steam cleaner.

6. Oxidation

- 10 min without vacuum at 800° C.

7. Porcelain application

- Porcelain applications are subject to the manufacturers instructions.
- PX GOLD SUPRA is compatible with low fusing ceramics (like Duceragold). Please observe the required coefficient of expansion.

8. Soldering

- Pre-bond: use PX Sold LF 950 solder at 950° C with a flux.
- Post-bond: use PX Sold 700 solder at 750° C with a flux.

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