N° 7350

**PX GOLD PLUS****PX DENTAL SA**

<b>JAUNE CLAIR</b>	<b>ALLIAGE DENTAIRE UNIVERSEL</b>
<b>HELL GELB</b>	<b>DENTAL UNIVERSALLEGIERUNGEN</b>
<b>GIALLO CHIARO</b>	<b>LEGA DENTALE UNIVERSAL</b>
<b>LIGHT YELLOW</b>	<b>DENTAL UNIVERSAL ALLOY</b>
<b>AMARILLO CLARO</b>	<b>ALEACIÓN DENTAL UNIVERSAL</b>

% = Au: 75.2, Pt: 9.0, Ag: 9.2, Cu: 3.0, Zn: 2.5, In: 1.0,  
Ir: < 1%

16.6 g/cm<sup>3</sup>

920-970° C

16.7/17.0 µm/m<sup>2</sup> C

ISO 9693 /22674

Fabriqu  par / Hergestellt von / Fabricato da /  
Made by / Fabricado por

**PX DENTAL SA**

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## 1. Waxing up, spruing, investing

- Follow the general prescriptions for bonding alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- Use only a high-heat phosphate-bonded investment material. For the alloy PX GOLD PLUS we recommend to use the PX Extravest investment.

## 2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 750° C for a soaking period of 30-45 min.

## 3. Crucibles and casting buttons

- May be cast in a graphite or ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

## 4. Melting

- For PX GOLD PLUS casting temperature is 1230° C with a graphite crucible, 1280° C with a ceramic crucible.
- The use of crystallized boric acid is recommended.
- Bench cool the casting ring.

## 5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and preceding oxidation clean the frameworks in an ultrasonic unit or with a steam cleaner.

## 6. Oxidation

- 10 min with vacuum at 800° C.

## 7. Porcelain application

- Porcelain applications are subject to the manufacturers instructions.
- PX GOLD PLUS is compatible with low fusing ceramics (like Duceragold). Please observe the required coefficient of expansion.

## 8. Soldering

- Pre-bond: use PX Sold LF 950 solder at 950° C with a flux.
- Post-bond: use PX Sold 700 solder at 750° C with a flux.

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