PXDENTAL SA

JAUNE	ALLIAGE DENTAIRE CÉRAMO MÉTALLIQUE
GELB	DENTAL AUFBRENNLEGIERUNGEN
GIALLO	LEGA DENTALE PER CERAMICA
YELLOW	DENTAL BONDING ALLOY
AMARILLO	ALECIÓN DENTAL CERAMO-METALLICA

% = Au: 86.5, Pt: 11.5, Zn: 1.4 ; Ir, Rh, Fe, Mn, Ta: < 1 %

18.9 g/cm³

1060-1140° C

14.3/14.6 µm/m° C

ISO 9693

Fabriqué par / Hergestellt von / Fabricato da /
Made by / Fabricado por

PX DENTAL SA

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1. Waxing up, spruing, investing

- Follow the general prescriptions for bonding alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- Use only a high-heat phosphate-bonded investment material. For the alloy PX PREMIUM BIO we recommend to use the PX Extravest investment.

2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 750-800° C for a soaking period of 30-45 min.

3. Crucibles and casting buttons

- Use a ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

4. Melting

- For PX PREMIUM BIO casting temperature is 1340° C.
- The use of crystallized boric acid is recommended.
- Bench cool the casting ring.

5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and preceding oxidation clean the frameworks in an ultrasonic unit or with a steam cleaner.

6. Oxidation

- 10 min under vacuum at 950° C.

7. Porcelain application

- Porcelain applications are subject to the manufacturers instructions.
- PX PREMIUM BIO is compatible with most known porcelain brands (high fusing).

8. Soldering

- Pre-bond: use PX Presold 1040 solder at 1040° C with a flux.
- Post-bond: use PX Postsold 800 solder at 850° C with a flux.

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