**BRASURE JAUNE APRES CUISSON
GELBES LOT NACH DEM BRAND
YELLOW SOLDER AFTER BONDING
SALDATURA GIALLA SECONDARIA**

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**Yellow solder after bonding
for precious ceramic alloys**

Composition in %:
Au: 75.0, Ag: 9.0, Cu: 7.0, Zn: 7.5, In: 1.5

Instructions for use

IMPORTANT DATA

Melting range 750 – 800°C
Working temperature 850°C

FOR THE FURNACE SOLDERING TECHNIQUE

- The surfaces are held together with wax, resin or composite. The gingival surface of the porcelain should be covered with wax 2 mm thick.
- The units should be parallel with a gap of 0.1-0.2 mm and a minimum soldering area of 6 mm².
- Make the soldering block as small as possible with a sufficient rigidity. All contact of the porcelain with investment should be avoided.
- Eliminate the wax with boiling water, the resin or composite with a Bunsen. Degrease the contact surfaces if necessary.
- Apply a flux to the cold surfaces as well as to the solder and place it between the solder joints.
- Dry the model for 6 min. at 400°C in the furnace with open door. Close the door and preheat the furnace at 500°C during at least 15 min. or 3 min. per unit.
- Switch on the vacuum and raise the temperature from 500°C to 850°C with a low heating rate of 55° C/min.
- Switch off the vacuum and maintain this temperature for 1 min. Remove the solder model and allow cooling depending the ceramic and the alloy

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