# PXDENTAL SA

BLANC	ALLIAGE DENTAIRE CÉRAMO MÉTALLIQUE
WEISS	DENTAL AUFBRENNLEGIERUNGEN
BIANCO	LEGA DENTALE PER CERAMICA
WHITE	DENTAL BONDING ALLOY
BLANCO	ALECION DENTAL CERAMO-METALLICA

% = Au: 51.2, Pd: 38.6, In: 8.6, Ga: 1.5, Ru: < 1%

14.5 g/cm<sup>3</sup>  
1180-1300° C  
13.8/14.1 µm/m<sup>2</sup> C

ISO 9693

Fabriqué par / Hergestellt von / Fabricato da /  
Made by / Fabricado por

**PX DENTAL SA**

Champs-Montants 16a, 2074 Marin  
T. +41 (0)32 924 21 20 F +41 (0)32 924 21 29  
www.pxdental.com

**PXDENTAL SA**  
CH - 2074 MARIN

## Mode d'emploi

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## 1. Waxing up, spruing, investing

- Follow the general prescriptions for casting alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- All regular or phosphate-bonded investments for precious metal alloys may be used. For the alloy PX BOND SF we recommend to use the PX Extravest investment.

## 2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 850° C for a soaking period of 30-45 min.

## 3. Crucibles and casting buttons

- PX BOND SF may be cast in a graphite, quartz or ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

## 4. Melting

- For PX BOND SF casting temperature is 1440° C with a ceramic crucible.
- The use of a pinch of boric acid is recommended.
- Bench cool the casting ring.

## 5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Finishing: use abrasive grinding points of your choice. Clean the frameworks in an ultrasonic unit or with a steam cleaner.

## 6. Oxidation

- 10 min without vacuum at 950°C.

## 7. Porcelain application

- Porcelain applications are subject to the manufacturer's instructions.
- PX BOND SF is compatible with most known porcelain brands. (high fusing).

## 8. Soldering

- Pre-bond : use PX PRESOLD 1110 at 1110° C with flux.
- Post-bond : use PX POSTSOLD 800 at 850° C with flux.

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