
ECCO SYSTEM Mode d'emploi Gebrauchsanweisung Istruzioni per l'uso

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
ECCO SYSTEM Instructions for Use Modo de empleo

Produit for Professional use only. Before begin please read "Important Notes on Use and Operation of Dental Attachments"

Description

Classic, proven snap-in anchor, low height, with 4 different male parts, for soldering, for connection casting and for casting.

Indications

Removable prostheses on root caps and implants, for example - total dentures

- intercalated prostheses and mixed intercalated and extended prostheses.

Contraindications

Unilateral free-end dentures without transverse connection

1. Positioning the male part

Using the parallelometer mandrin, position the male part on the flat milled root cap (note the position key) and fix it with adhesive wax. Design the soldering block as small as possible and ensure easy access to the flame.

2. Hardening of the male part

The Au-Ag male part is the most abrasion resistant. It must be annealed (400°C / 15 min.) to reach a hardness of 345 Vickers.

The Au-Pt male parts are self hardening. The male parts welded on the root caps are to be cooled slowly to room temperature.

HMA is the "softest" alloy for the male parts. They can be out-casted. Use HMA males only if no other solution is possible. After casting, cool slowly to room temperature.

When using PMMA burnout resin males, use an alloy with a 0.2% yield strength of at least 500 N/mm².

3. Polishing

When polishing, be careful not to damage the male part. Protect it with the auxiliary part.

4. Assembling the female part

Assemble with the inner mounting sockets (A) Screw the inner mounting sockets with a torque wrench into the retaining cap. Place the complete female parts onto the male parts. The female parts must be placed parallel and axially in the direction of insertion. Lock the female parts and the undercuts. Then coat and cure the prosthesis. Remove the excess resin from the bottom of the female parts. Unscrew the inner mounting sockets and replace them with the original inner sockets. Assembly with the original inner females

Lightly coat the original female parts with Vaseline. Cover the original female parts (Au-Pt) with a plastic ring. Place the original female parts on the male parts. The female parts must be parallel and axially in the direction of insertion. Fix the female parts and the undercuts. Then coat and cure the prosthesis.

5. Assembling the female part in the mouth

Place the auxiliary parts absolutely parallel and axially in the direction of insertion on the model. Polymerize the prosthesis with the auxiliary parts. Then remove them. The resulting hole should be slightly larger than the original female part.

In the mouth: place a perforated dam over the gingiva and the male parts. The female parts can be assembled (assembly or original) as described in section 4. Block out the undercuts with wax or silicone. Place a small amount of self-curing resin in the cavities of the prosthesis and position the prosthesis in situ. After the resin has cured, remove the prosthesis, add more resin or remove the excess.

6. Removing the female part

Unscrew the inner female part. Warm up the extractor and screw it into the retention cap (warm up the resin). Use the extractor to remove the resin retention cap from the denture.

7. Check-up

Caution: Resistant resin elements in attachments have enormous advantages.

However, the elasticity of all resins leads to faster wear in the oral cavity. In addition, abrasive substances (e.g. toothpaste) can damage the resin. For these reasons, maximum hygiene is required and all resin elements should be replaced at least once a year.

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