N° 7301

PX Bio Gold LF

**PX Dental SA**Bvd Eplatures 42, c/o PX group  
CH-2304 La Chaux-de-Fonds**PX DENTAL**

<b>JAUNE</b>	<b>ALLIAGE DENTAIRE UNIVERSEL</b>
<b>GELB</b>	<b>DENTAL UNIVERSALLEGIERUNGEN</b>
<b>GIALLO</b>	<b>LEGA DENTALE UNIVERSAL</b>
<b>YELLOW</b>	<b>DENTAL UNIVERSAL ALLOY</b>
<b>AMARILLO</b>	<b>ALEACIÓN DENTAL UNIVERSAL</b>

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## Instructions for use

Before you start, read the "Important instructions for use and processing of the Dental Alloys" Alloys are products for professional use only.

### Indications

Central or multisurface inlays, ¾ crowns, full crowns, Short or middle span bridges

### Contraindications

Do not use in know hypersensitivity to one or several metals contained in the alloy, Avoid occlusal and proximal contacts of different alloys, Unfavorable periodontal environment

### 1. Waxing up, spruing, investing

- Follow the general prescriptions for bonding alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- Use only a high-heat phosphate-bonded investment material.

### 2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 700° C for a soaking period of 20-45 min.

### 3. Crucibles and casting buttons

- May be cast in a graphite or ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

### 4. Melting

- For PX BIO GOLD LF casting temperature is 1210° C with a graphite crucible, 1260° C with a ceramic crucible.
- The use of crystallized boric acid is recommended.
- Bench cool the casting ring.

### 5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and preceding oxidation clean the frameworks in an ultrasonic unit or with a steam cleaner.

### 6. Oxidation

- 10 min without vacuum at 800° C.

### 7. Porcelain application

- Porcelain applications are subject to the manufacturers instructions.
- PX BIO GOLD LF is compatible with low fusing ceramics (like Duceragold). Please observe the required coefficient of expansion.

### 8.Soldering

- Pre-bond: use PX Sold LF 950 solder at 950° C with a flux.
- Post-bond: use PX Sold 700 solder at 750° C with a flux.

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