N° 7120

PX Bond SF

**PX Dental SA**

 Bvd Eplatures 42, c/o PX group  
 CH-2304 La Chaux-de-Fonds

**PX DENTAL**

BLANC ALLIAGE DENTAIRE CÉRAMO MÉTALLIQUE  
 WEISS DENTAL AUFBRENNLEGIERUNGEN  
 BIANCO LEGA DENTALE PER CERAMICA  
 WHITE DENTAL BONDING ALLOY  
 BLANCO ALEACION DENTAL CERAMO-METALICA

**Business address:**

 Champs-Montants 16a  
 CH-2074 Marin

 T. +41(0)32 924 21 20  
 pxdental@pxgroup.com

% = Au: 51.2, Pd: 38.6, In: 8.6, Ga: 1.5, Ru: &lt; 1%

 14.5 g/cm<sup>3</sup>

1180-1300°C

13.8/14.1 µm/m°C

Type 4 selon ISO 22674

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## Instructions for use

### Indications

¾ crowns, full crowns, Short and long span bridges, milled telescope crowns, implant suprastructures.

### Contraindications

Do not use in know hypersensitivity to one or several metals contained in the alloy, Avoid occlusal and proximal contacts of different alloys, Unfavorable periodontal environment

#### 1. Waxing up, spruing, investing

- Follow the general prescriptions for casting alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- All regular or phosphate-bonded investments for precious metal alloys may be used. For the alloy PX BOND SF we recommend to use the PX Extravest investment.

#### 2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.

- Keep the cylinder at 850° C for a soaking period of 30-45 min.

#### 3. Crucibles and casting buttons

- PX BOND SF may be cast in a graphite, quartz or ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

#### 4. Melting

- For PX BOND SF casting temperature is 1440° C with a ceramic crucible.
- The use of a pinch of boric acid is recommended.
- Bench cool the casting ring.

#### 5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Finishing: use abrasive grinding points of your choice. Clean the frameworks in an ultrasonic unit or with a steam cleaner.

#### 6. Oxidation

- 10 min without vacuum at 950°C.

#### 7. Porcelain application

- Porcelain applications are subject to the manufacturer's instructions.
- PX BOND SF is compatible with most known porcelain brands. (high fusing).

#### 8. Soldering

- Pre-bond : use PX PRESOLD 1110 at 1110° C with flux.
- Post-bond : use PX POSTSOLD 800 at 850° C with flux.

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