N° 7133

PX Bondaloy

PX Dental SA

 Bvd Eplatures 42, c/o PX group
 CH-2304 La Chaux-de-Fonds

PX DENTAL

BLANC ALLIAGE DENTAIRE CÉRAMO MÉTALLIQUE
 WEISS DENTAL AUFBRENNLEGIERUNGEN
 BIANCO LEGA DENTALE PER CERAMICA
 WHITE DENTAL BONDING ALLOY
 BLANCO ALEACION DENTAL CERAMO-METALICA

Business address:

 Champs-Montants 16a
 CH-2074 Marin

 T. +41(0)32 924 21 20
 pxdental@pxgroup.com

% = Au: 54.2, Pd: 31.0, Ag: 4.8, In: 9.0, Ru, Ga: < 1%

 14.6 g/cm³
 1160-1260° C
 14.3/14.5 µm/m° C

Type 4 selon ISO 22674

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Instructions for use

Indications

¾ crowns, full crowns, Short and long span bridges, milled telescope crowns, implant suprastructures.

Contraindications

Do not use in know hypersensitivity to one or several metals contained in the alloy, Avoid occlusal and proximal contacts of different alloys, Unfavorable periodontal environment

1. Waxing up, spruing, investing

- Follow the general prescriptions for casting alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- All regular or phosphate-bonded investments for precious metal alloys may be used. For the alloy PX BONDALLOY we recommend to use the PX Extravest investment.

2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.

- Keep the cylinder at 850° C for a soaking period of 30-45 min.

3. Crucibles and casting buttons

- PX BONDALLOY may be cast in a graphite, quartz or ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

4. Melting

- For PX BONDALLOY casting temperature is 1440° C with a ceramic crucible.
- The use of a pinch of boric acid is recommended.
- Bench cool the casting ring.

5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Finishing: use abrasive grinding points of your choice. Clean the frameworks in an ultrasonic unit or with a steam cleaner.

6. Oxidation

- 10 min without vacuum at 950°C.

7. Porcelain application

- Porcelain applications are subject to the manufacturer's instructions.
- PX BONDALLOY is compatible with most known porcelain brands. (high fusing).

8. Soldering

- Pre-bond : use PX PRESOLD 1110 at 1110° C with flux.
- Post-bond : use PX POSTSOLD 800 at 850° C with flux.

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