N° 7171

PX Ecopal

**PX Dental SA**

Bvd Eplatures 42, c/o PX group  
CH-2304 La Chaux-de-Fonds

**PX DENTAL**

BLANC	ALLIAGE DENTAIRE CÉRAMO MÉTALLIQUE
WEISS	DENTAL AUFBRENNLEGIERUNGEN
BIANCO	LEGA DENTALE PER CERAMICA
WHITE	DENTAL BONDING ALLOY
BLANCO	ALEACION DENTAL CERAMO-METALICA

**Business address:**

Champs-Montants 16a  
CH-2074 Marin

T. +41(0)32 924 21 20  
pxdental@pxgroup.com

% = Pd: 61.4, Ag: 24.5, Zn: 2.0, In: 2.0, Sn: 10.0,  
Ru, B: < 1%

11.3 g/cm<sup>3</sup>  
1180-1230° C  
14.5/14.9 µm/m°C

Type 4 selon ISO 22674

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## Instructions for use

### Indications

¾ crowns, full crowns, Short and long span bridges, milled telescope crowns, implant suprastructures, Cast partial dental frames

### Contraindications

Do not use in know hypersensitivity to one or several metals contained in the alloy, Avoid occlusal and proximal contacts of different alloys, Unfavorable periodontal environment

#### 1. Waxing up, spruing, investing

- Follow the general prescriptions for bonding alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- Use only a high-heat phosphate-bonded investment material.

#### 2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 850° C for a soaking period of 20-45 min.

#### 3. Crucibles and casting buttons

- Use a ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

#### 4. Melting

- For PX ECOPAL casting temperature is 1430° C.
- The use of crystallized boric acid is recommended.
- Bench cool the casting ring.

#### 5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and preceding oxidation clean the frameworks in an ultrasonic unit or with a steam cleaner.

#### 6. Oxidation

- 10 min without vacuum at 950° C.

#### 7. Porcelain application

- Porcelain applications are subject to the manufacturers instructions.
- PX ECOPAL is compatible with most known porcelain brands (high fusing).

#### 8. Soldering

- Pre-bond: use PX PRESOLD 1110 solder at 1110° C with a flux.
- Post-bond: use PX POSTSOLD 800 solder at 850° C with a flux.

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