N° 7303

PX Gold Plus

**PX Dental SA**Bvd Eplatures 42, c/o PX group
CH-2304 La Chaux-de-Fonds**PX DENTAL**

JAUNE CLAIR	ALLIAGE DENTAIRE UNIVERSEL
HELL GELB	DENTAL UNIVERSALLEGIERUNGEN
GIALLO CHIARO	LEGA DENTALE UNIVERSAL
LIGHT YELLOW	DENTAL UNIVERSAL ALLOY
AMARILLO CLARO	ALEACIÓN DENTAL UNIVERSAL

Business address:Champs-Montants 16a
CH-2074 MarinT. +41(0)32 924 21 20
pxdental@pxgroup.com% = Au: 75.2, Pt: 9.0, Ag: 9.2, Cu: 3.0, Zn: 2.5, In: 1.0,
Ir: < 1%16.6 g/cm³
920-970°C
16.7/17.0 µm/m°C

Type 4 selon ISO 22674

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Instructions for use

Indications

Central or multisurface inlays, ¾ crowns, full crowns, Short and long span bridges, milled telescope crowns, implant suprastructures.

Contraindications

Do not use in know hypersensitivity to one or several metals contained in the alloy, Avoid occlusal and proximal contacts of different alloys, Unfavorable periodontal environment

1. Waxing up, spruing, investing

- Follow the general prescriptions for bonding alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- Use only a high-heat phosphate-bonded investment material.

For the alloy PX GOLD PLUS we recommend to use the PX Extravest investment.

2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 750° C for a soaking period of 30-45 min.

3. Crucibles and casting buttons

- May be cast in a graphite or ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

4. Melting

- For PX GOLD PLUS casting temperature is 1230° C with a graphite crucible, 1280° C with a ceramic crucible.

- The use of crystallized boric acid is recommended.
- Bench cool the casting ring.

5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and preceding oxidation clean the frameworks in an ultrasonic unit or with a steam cleaner.

6. Oxidation

- 10 min with vacuum at 800° C.

7. Porcelain application

- Porcelain applications are subject to the manufacturers instructions.
- PX GOLD PLUS is compatible with low fusing ceramics (like Duceragold). Please observe the required coefficient of expansion.

8. Soldering

- Pre-bond: use PX Sold LF 950 solder at 950° C with a flux.
- Post-bond: use PX Sold 700 solder at 750° C with a flux.

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