N° 7220

PX Palfit

PX Dental SA

Bvd Eplatures 42, c/o PX group
CH-2304 La Chaux-de-Fonds

PX DENTAL

BLANC	ALLIAGE DENTAIRE CÉRAMO MÉTALLIQUE
WEISS	DENTAL AUFBRENNLEGIERUNGEN
BIANCO	LEGA DENTALE PER CERAMICA
WHITE	DENTAL BONDING ALLOY
BLANCO	ALEACION DENTAL CERAMO-METALICA

Business address:

Champs-Montants 16a
CH-2074 Marin

T. +41(0)32 924 21 20
pxdental@pxgroup.com

% = Pd: 53.8, Ag: 38.0, Sn: 8.0, Ru: < 1%

11.4 g/cm³

1190-1270°C

14.6/14.9 µm/m°C

Type 4 selon ISO 22674

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Instructions for use

Indications

¾ crowns, full crowns, Short and long span bridges, milled telescope crowns, implant suprastructures, Cast partial dental frames

Contraindications

Do not use in know hypersensitivity to one or several metals contained in the alloy, Avoid occlusal and proximal contacts of different alloys, Unfavorable periodontal environment

1. Waxing up, spruing, investing

- Follow the general prescriptions for bonding alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- Use only a high-heat phosphate-bonded investment material.

2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 850° C for a soaking period of 20-45 min.

3. Crucibles and casting buttons

- Use a ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

4. Melting

- For PX PALFIT casting temperature is 1420° C.
- The use of crystallized boric acid is recommended.
- Bench cool the casting ring.

5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and preceding oxidation clean the frameworks in an ultrasonic unit or with a steam cleaner.

6. Oxidation

- 10 min without vacuum at 950° C.

7. Porcelain application

- Porcelain applications are subject to the manufacturers instructions.
- PX PALFIT is compatible with most known porcelain brands (high fusing).

8. Soldering

- Pre-bond: use PX PRESOLD 1110 solder at 1110° C with a flux.
- Post-bond: use PX POSTSOLD 800 solder at 850° C with a flux.

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