N° 7040

PX Premium C

**PX Dental SA**

Bvd Eplatures 42, c/o PX group  
CH-2304 La Chaux-de-Fonds

**PX DENTAL**

BLANC	ALLIAGE DENTAIRE CÉRAMO MÉTALLIQUE
WEISS	DENTAL AUFBRENNLEGIERUNGEN
BIANCO	LEGA DENTALE PER CERAMICA
WHITE	DENTAL BONDING ALLOY
BLANCO	ALEACION DENTAL CERAMO-METALICA

**Business address:**

Champs-Montants 16a  
CH-2074 Marin

T. +41(0)32 924 21 20  
pxdental@pxgroup.com

% = Au: 75.0; Pd: 18.5; Ag: 1.5; In: 2.0; Sn: 2.0, Ir, Ru,  
Cu, Zn : < 1%

16.1 g/cm<sup>3</sup>

1120 - 1240°C

14.1 / 14.4 μm/m°C

Type 4 selon ISO 22674

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## Instructions for use

### Indications

Central or multisurface inlays, ¾ crowns, full crowns, Short and long span bridges, milled telescope crowns, implant suprastructures.

### Contraindications

Do not use in know hypersensitivity to one or several metals contained in the alloy, Avoid occlusal and proximal contacts of different alloys, Unfavorable periodontal environment

1. Waxing up, spruing, investing
  - Follow the general prescriptions for bonding alloys.
  - The ultimate thickness of the cast coping should not be less than 0.3 mm.

• Use only a high-heat phosphate-bonded investment material.

For the alloy PX PREMIUM C we recommend to use the PX Extravest investment.

### 2. Preheating

• Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.

• Keep the cylinder at 850° C for a soaking period of 30-45 min.

### 3. Crucibles and casting buttons

• The use of a ceramic crucible is recommended.

• Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

### 4. Melting

• For PX PREMIUM C casting temperature is 1390° C with a graphite crucible, 1440° C with a ceramic crucible.

• The use of crystallized boric acid is recommended.

• Bench cool the casting ring.

### 5. Treatment of the casting

• Devesting: by blasting (grain size of 50-125 µm).

• Following casting or soldering, the frame may be pickled.

• Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and preceding oxidation clean the frameworks in an ultrasonic unit or with a steam cleaner.

### 6. Oxidation

• 10 min under vacuum at 950° C.

### 7. Porcelain application

• Porcelain applications are subject to the manufacturers instructions.

• PX PREMIUM C is compatible with most known porcelain brands (high fusing).

### 8. Soldering

• Pre-bond: use PX PRESOLD 1100 solder at 1100° C with a flux.

• Post-bond: use PX POSTSOLD 800 solder at 850° C with a flux.

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