N° 7030

PX Premium S

**PX Dental SA**Bvd Eplatures 42, c/o PX group  
CH-2304 La Chaux-de-Fonds**PX DENTAL**

JAUNE CLAIR	ALLIAGE DENTAIRE CÉRAMO MÉTALLIQUE
HELLGELB	DENTAL AUFBRENNLEGIERUNGEN
GIALLO CHIARO	LEGA DENTALE PER CERAMICA
LIGHT YELLOW	DENTAL BONDING ALLOY
AMARILLO CLARO	ALEACION DENTAL CERAMO-METALICA

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## Instructions for use

### Indications

Central or multisurface inlays, ¾ crowns, full crowns, Short and long span bridges, milled telescope crowns, implant suprastructures.

### Contraindications

Do not use in know hypersensitivity to one or several metals contained in the alloy, Avoid occlusal and proximal contacts of different alloys, Unfavorable periodontal environment

#### 1. Waxing up, spruing, investing

- Follow the general prescriptions for bonding alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- Use only a high-heat phosphate-bonded investment material.

For the alloy PX PREMIUM S we recommend to use the PX Extravest investment.

#### 2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.

- Keep the cylinder at 850° C for a soaking period of 30-45 min.

#### 3. Crucibles and casting buttons

- The use of a ceramic crucible is recommended.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

#### 4. Melting

- For PX PREMIUM S casting temperature is 1340° C with a graphite crucible, 1390° C with a ceramic crucible
- The use of crystallized boric acid is recommended.
- Bench cool the casting ring.

#### 5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and preceding oxidation clean the frameworks in an ultrasonic unit or with a steam cleaner.

#### 6. Oxidation

- 10 min under vacuum at 950° C.

#### 7. Porcelain application

- Porcelain applications are subject to the manufacturers instructions.
- PX PREMIUM S is compatible with most known porcelain brands (high fusing).

#### 8. Soldering

- Pre-bond: use PX PRESOLD 1100 solder at 1100° C with a flux.
- Post-bond: use PX POSTSOLD 800 solder at 850° C with a flux.

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