N° 7105

PX Premium Yellow

PX Dental SA

 Bvd Eplatures 42, c/o PX group
 CH-2304 La Chaux-de-Fonds

PX DENTAL

JAUNE ALLIAGE DENTAIRE CÉRAMO MÉTALLIQUE
 GELB DENTAL AUFBRENNLEGIERUNGEN
 GIALLO LEGA DENTALE PER CERAMICA
 YELLOW DENTAL BONDING ALLOY
 AMARILLO ALEACION DENTAL CERAMO-METALICA

Business address:

 Champs-Montants 16a
 CH-2074 Marin

 T. +41(0)32 924 21 20
 pxdental@pxgroup.com

% = Au: 76.5, Pt: 13.5, Ag: 6.5, Zn: 2.5, Ta, Ir, Ru: < 1%

17.8 g/cm³

1020-1120°C

14.4/14.8 µm/m°C

Type 4 selon ISO 22674

Mode d'emploi

Indications :

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Instructions for use

Indications

Central or multisurface inlays, ¾ crowns, full crowns, Short and long span bridges, milled telescope crowns, implant suprastructures.

Contraindications

Do not use in know hypersensitivity to one or several metals contained in the alloy, Avoid occlusal and proximal contacts of different alloys, Unfavorable periodontal environment

1. Waxing up, spruing, investing

- Follow the general prescriptions for bonding alloys.
- The ultimate thickness of the cast coping should not be less than 0.3 mm.
- Use only a high-heat phosphate-bonded investment material.

For the alloy PX PREMIUM YELLOW we recommend to use the PX Extravest investment.

2. Preheating

- Observe manufacturer's recommendations with regard to setting times, temperature levels, etc.
- Keep the cylinder at 750-800° C for a soaking period of 30-45 min.

3. Crucibles and casting buttons

- Use a ceramic crucible.
- Add at least 1/3 of new alloy. When remelted casting buttons should be perfectly clean.

4. Melting

- For PX PREMIUM YELLOW casting temperature is 1320° C.
- The use of crystallized boric acid is recommended.
- Bench cool the casting ring.

5. Treatment of the casting

- Devesting: by blasting (grain size of 50-125 µm).
- Following casting or soldering, the frame may be pickled.
- Preparing the framework and cleaning: trim the framework preferably using carbide burs or fine carborundum stones at low speed another blasting (grain size of 50-125 µm) and preceding oxidation clean the frameworks in an ultrasonic unit or with a steam cleaner.

6. Oxidation

- 10 min under vacuum at 900° C.

7. Porcelain application

- Porcelain applications are subject to the manufacturers instructions.
- PX PREMIUM YELLOW is compatible with most known porcelain brands (high fusing).

8. Soldering

- Pre-bond: use PX SOLD LF 950 solder at 950° C with a flux.
- Post-bond: use PX SOLD 740 solder at 790° C with a flux.

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