
MINI-NPG

Mode d'emploi

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

MINI-NPG

Instructions for Use Modo de empleo

Produit for Professional use only. Before begin please read "Important Notes on Use and Operation of Dental Attachments"

Description

Extracoronar friction slide, consisting of a male part made of precious metal alloy for casting fittings or of burnable resin and female parts made of steel. Due to the stable construction a milled contouring is not necessary.

Indications

Rigid partial dentures with dento-mucosal support

Contraindications

Prostheses with unilateral extension saddle, without transverse connection

1. Cleaning

Clean the female and male parts, which must be free of grease and other impurities (e.g. with a steam jet).

2. Positioning of the male part

Place the male part in the parallelogram rod and insert it into the wax model.

Avoid any contact between the ceramic and the HMA male part suitable for the fitting casting (incompatible TEC)

3. Groove

A fine groove between the male part and the wax prevents possible overflow of the casting alloy onto the slide. Do not apply any surfactant to the slide elements.

When using the male part, make sure to use an alloy with a 0.2% yield strength of at least 500 MPa.

4. Preheating temperatures

Conventional alloys 700°C

Ceramic-metal alloys MP 900°C

5. Pickling

After casting, remove from mould and pickle in hot 10% sulphuric acid or acid substitute. Do not sandblast.

6. Slide attachment fit

Check the attachments for a smooth fit.

7. Mounting in resin

On the model, insert the female part on the male part, fill all the openings and all the retentive parts. Carry out the assembly and proceed with the polymerization.

8. Bonding

The female part can also be bonded into the stellite framework

9. Relining

We only recommend the indirect method, which consists of making an impression and fabricating a model: to do this, insert the transfer pin into the impression, fabricate the model and finish the work.

10. Disassembling and assembling the sleeve

Remove the sleeve with an instrument (scalpel). Place the sleeve on the mounting instrument and press it into the empty female part until it snaps into place perceptibly.

11. Warning concerning resin elements

Strong resin elements in attachments have enormous advantages. However, the elasticity of all resins leads to faster wear in the oral cavity. Therefore, resin elements should be replaced at least once a year.

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