## PX Dental SA

Bvd Eplatures 42, c/o PX group  
CH-2304 La Chaux-de-Fonds  
www.pxdental.comUtilisation pour les alliages céramo-métalliques / Geeignet für Edelmetall  
Aufbrennlegierungen  
For precious ceramic alloys / per le leghe per ceramica**Business address:**  
Champs-Montants 16a  
CH-2024 Marin  
T. +41(0)32 924 21 20  
pxdental@pxgroup.com

## BRASURE JAUNE APRES CUISSON

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## Instructions for use

## Important data

Melting range: 750 – 800°C  
Working temperature: 850°C*Before you start, read the "Important instructions for use and processing of the Dental Alloys". Alloys are products for professional use only.***Contraindications***Do not use in know hypersensitivity to one or several metals contained in the alloy, Avoid occlusal and proximal contacts of different alloys, Unfavorable periodontal environment*

## FOR THE FURNACE SOLDERING TECHNIQUE

- The surfaces are held together with wax, resin or composite. The gingival surface of the porcelain should be covered with wax 2 mm thick.
- The units should be parallel with a gap of 0.1-0.2 mm and a minimum soldering area of 6 mm<sup>2</sup>.
- Make the soldering block as small as possible with a sufficient rigidity. All contact of the porcelain with investment should be avoided.
- Eliminate the wax with boiling water, the resin or composite with a Bunsen. Degrease the contact surfaces if necessary.
- Apply a flux to the cold surfaces as well as to the solder and place it between the solder joints.
- Dry the model for 6 min. at 400°C in the furnace with open door. Close the door and preheat the furnace at 500°C during at least 15 min. or 3 min. per unit.
- Switch on the vacuum and raise the temperature from 500°C to 850°C with a low heating rate of 55° C/min.

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