## PX Dental SA

Bvd Eplatures 42, c/o PX group  
CH-2304 La Chaux-de-Fonds  
www.pxdental.com

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## LIGHT-YELLOW SOLDER BEFORE BONDING

For the alloys : :  
PX Bond, SF, Special, Budget, Extrapall, Ecopal, PalFit

Composition in %  
Au: 85.0, Pd: 10.5, Zn: 4.5

## Instructions for use

## Important data

Melting range: 1030 – 1110°C  
Working temperature: 1110°C

*Before you start, read the "Important instructions for use and processing of the Dental Alloys"  
Alloys are products for professional use only.*

*Contraindications*

*Do not use in know hypersensitivity to one or several metals contained in the alloy, Avoid occlusal and proximal contacts of different alloys, Unfavorable periodontal environment*

- The surfaces are held together with wax, resin or composite.
- The units should be parallel with a gap of 0.1-0.2 mm and a minimum soldering area of 6 mm<sup>2</sup>.
- Make the soldering block as small as possible with a sufficient rigidity.
- Eliminate the wax with boiling water, the resin or composite with an open flame. Degrease the contact surfaces if necessary.
- Apply a flux to the cold surfaces as well as to the solder.
- Heat the solder block and the surfaces to be soldered to working temperature (1110°C).
- Apply the solder only when the working temperature is reached.
- Bench cool the solder and remove the soldering block.
- Pickle with a usual pickling liquid.

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