PX Dental SA

Bvd Eplatures 42, c/o PX group
CH-2304 La Chaux-de-Fonds
www.pxdental.com

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Solder for combined works between a gold based alloy and a NI-CR or CO-CR alloy

Composition in %
Pd: 30.0, Ag: 46.0, Cu: 24.0

Instructions for use

Important data

Melting range: 860 – 940°C
Working temperature: 940°C

Before you start, read the "Important instructions for use and processing of the Dental Alloys" Alloys are products for professional use only.

Contraindications:

Do not use in know hypersensitivity to one or several metals contained in the alloy, Avoid occlusal and proximal contacts of different alloys, Unfavorable periodontal environment Attention: this product contains nickel which may have allergizing effects or provoke localised reactions. Please consider the use of the product and the particular sensitivity of the patient.

- The surfaces are held together with wax, resin or composite.
- The units should be parallel with a gap of 0.1-0.2 mm and a minimum soldering area of 6 mm².
- Make the soldering block as small as possible with a sufficient rigidity.
- Eliminate the wax with boiling water, the resin or composite with an open flame. Degrease the contact surfaces if necessary.
- Apply a flux to the cold surfaces as well as to the solder.
- Heat the solder block and the surfaces to be soldered to working temperature (940°C).
- Apply the solder only when the working temperature is reached.
- Bench cool the solder and remove the soldering block.
- Pickle with a usual pickling liquid.

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