
SFERA SYSTEM Mode d'emploi Gebruichsanweisung Istruzioni per l'uso

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

SFERA SYSTEM Instructions for Use Modo de empleo

Produit for Professional use only. Before begin please read "Important Notes on Use and Operation of Dental Attachments"

Description

Classic, proven snap-in anchor, low height, with 4 different male parts, for soldering, for connection casting and for casting.

Indications

Removable prostheses on root caps and implants, for example

- total dentures
- intercalated prostheses and mixed intercalated and extended prostheses.

Contraindications

Unilateral free-end dentures without transverse connection

1. Positioning the male part

Using the parallelometer mandrin, position the male part on the flat milled root cap (note the position key) and fix it with adhesive wax. Design the soldering block as small as possible and ensure easy access to the flame.

2. Hardening of the male part

The Au-Pt male parts are self hardening. The male parts welded on the root caps are to be cooled slowly to room temperature. HMA can be out-casted. After casting, cool slowly to room temperature. When using PMMA burnout resin males, use an alloy with a 0.2% yield strength of at least 500 Mpa.

3. Polishing

When polishing, be careful not to damage the male part. Protect it with the auxiliary part.

4. Mounting the female part

Before curing, fill the retentive parts with e.g. silicone. The resin must not penetrate between the female and male parts. If there are several anchors, place the female parts parallel to each other. For assembly, position the spacer to allow the resilient solution between the male part between the male part, the root cap and the female part. Remove the spacer disc after the prosthesis has been fabricated. In the mouth, the prosthesis descends under load into the space created by the spacer.

5. Mounting the female part in the mouth

Complete the fabrication of the prosthesis with the auxiliary part and then remove it. The resulting cavity should be slightly larger than the original female part. Drill the dam and place it over the male part, the cap and the gingiva. Place the original female part over the male part. Fill the female part and insert the prosthesis. After the resin has cured, remove the prosthesis. If necessary, fill the remaining cavity with fresh resin and cure in a pressure pot if possible. Remove the excess resin first.

6. Temporary solution with the dummy

The auxiliary part can replace the original female part and serve as a temporary solution.

If there are several anchors, the force required to remove the prosthesis can be minimized. Before mounting in the resin, lightly grind the outer surfaces of the auxiliary part.

7. Activating the female part

Place the prosthesis on a stable support. Direct the activator axially into the female part by applying pressure and rotation to the lamellae to tighten them.

8. Deactivating the female part

Place the prosthesis on a stable surface. Direct the deactivator axially into the female part and press the plates apart.

9. Relining

Root caps: Place the transfer jig in the impression or in the female part of the prosthesis and fabricate the model. Implant work: Instead of using the Sfera transfer pin, use a replica or the analog appropriate for the implant.

10. Removing the female part

Insert the extractor into the female part. Heat the extractor, hold it with tweezers and pull the female part from the denture resin.

11. Mounting the plastic ring

Place the plastic rings on the stamp. Integrate the handle. Insert the mounting device into the female part. Push the handle until the ring is flush with the female part.

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