**PX Dental SA**

Bvd Eplatures 42, c/o PX group
CH-2304 La Chaux-de-Fonds
www.pxdental.com

Business address:

Champs-Montants 16a
CH-2024 Marin
T. +41(0)32 924 21 20
pxdental@pxgroup.com

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For the alloys : :
PX Gold, PX Inlaygold, PX Ecocast

Composition in %
Au: 75, Pd: 2.8, Ag: 1.5, Cu: 12.5, Zn: 2.9, In: 5.3

Instructions for use**Important data**

Melting range: 790 – 830°C
Working temperature: 830°C

*Before you start, read the "Important instructions for use and processing of the Dental Alloys"
Alloys are products for professional use only.*

Contraindications

Do not use in know hypersensitivity to one or several metals contained in the alloy, Avoid occlusal and proximal contacts of different alloys, Unfavorable periodontal environment

- The surfaces are held together with wax, resin or composite.
- The units should be parallel with a gap of 0.1-0.2 mm and a minimum soldering area of 6 mm².
- Make the soldering block as small as possible with a sufficient rigidity.
- Eliminate the wax with boiling water, the resin or composite with an open flame. Degrease the contact surfaces if necessary.
- Apply a flux to the cold surfaces as well as to the solder.
- Heat the solder block and the surfaces to be soldered to working temperature (830°C).
- Apply the solder only when the working temperature is reached.
- Bench cool the solder and remove the soldering block.
- Pickle with a usual pickling liquid.

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